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

**SUMMARY**

1 - THE INSTALLATION AND COMPONENTS OF EQUIPMENT.....	4
1.1 - INSTRUCTION.....	4
1.2 - COMPONENTS.....	5
1.2.1 - The following parts should be included in a complete UDS-E ultrasonic piezo scaler:	5
1.2.2 - Product performance and structure.....	6
1.2.3 - Scope of application	6
1.3 - THE MAIN TECHNICAL SPECIFICATIONS.....	6
1.4 - INSTALLATION OF THE MAIN COMPONENTS.....	7
1.4.1 - Sketch map for installation and connection.....	7
2 - PRODUCT FUNCTION AND USAGE.....	10
2.1 - SCALING FUNCTION	10
2.1.1 - Usage	10
2.1.2 - Instruction for main components of detachable handpiece (See Figure 6)	12
2.1.3 Instruction of torque wrench (See Figure 7)	12
2.2 - AUTO-WATER SUPPLY FUNCTION	13
2.3 - ENDO FUNCTION	14
3 - STERILIZATION AND MAINTENANCE.....	15
3.1 - STERILIZATION OF DETACHABLE HANDPIECE.....	15
3.1.1 - Autoclaved in high temperature/pressure:.....	15
3.1.2 - Notice	15
3.2 - STERILIZATION OF SCALING TIPS AND ENDOCHUCK.....	16
3.3 - STERILIZATION OF TORQUE WRENCH AND ENDO WRENCH	16
3.4 - CLEANING OF TIPS, ENDOCUCK, TORQUE WRENCH AND ENDO WRENCH	16
3.5 - TROUBLESHOOTING AND NOTES	17
3.5.1 - Troubleshooting.....	17
3.5.2 - Notes	19
4 - PRECAUTION.....	21
4.1 - NOTICE WHEN USING EQUIPMENT.....	21
4.2 - CONTRAINDICATION	22
4.3 - STORAGE AND MAINTENANCE	22
4.4 - TRANSPORTATION	22

4.5 - WORKING CONDITION	22
5 - AFTER SERVICE	23
6 - SYMBOLS INSTRUCTION	23
7 - ENVIRONMENTAL PROTECTION.....	24
8 - DECLARATION OF CONFORMITY	24
8.1 - PRODUCT CONFORMITY THE FOLLOWING STANDARDS:.....	24
9 - FOR TECHNICAL DATA, PLEASE CONTACT	24
9.1 - EMC - DECLARATION OF CONFORMITY	25

1 - THE INSTALLATION AND COMPONENTS OF EQUIPMENT

1.1 - INSTRUCTION

Carlo De Giorgi S.r.l. is a professional manufacturer in researching, developing and producing ultrasonic piezo scalers. The product is mainly used for teeth cleaning and also an indispensable equipment for tooth disease prevention and treatment. The new product UDS-E ultrasonic piezo scaler has scaling, perio, endo and auto-water supply functions, it contains the following features:

1. With auto-water supply function, easier operation.
2. The inner water pipe is made from antiseptic material, clinical fluids, such as Hydrogen peroxide, chlorhexidine and sodium Hypochlorite, etc., may be used perfectly under the mode of auto-water supply system so that the performance of perio and endo is improved dramatically.
3. The handpiece is detachable and can be autoclaved under the high temperature of 135°C and the pressure of 0.22MPa.
4. Automatic frequency tracking ensures that the machine always works on the best frequency and performs more steadily.
5. Controlled by PC, easy operation and more efficient for scaling.

These features make UDS-E become a new generation product in the world dental market.

1.2 - COMPONENTS

1.2.1 - The following parts should be included in a complete UDS-E ultrasonic piezo scaler:

Number	Description	Type
1	Main unit	234 mm × 193 mm × 109 mm
2	Detachable handpiece	Φ 18 mm × 117 mm
3	Adapter of power supply	133 mm × 58 mm × 35mm
4	Foot switch	Φ 92 mm × 26mm
5	Scaling tip	
6	Water pipe	Φ 6 mm × 4mm
7	Torque wrench	Φ 32 mm × 39 mm
8	Waterproof "o" ring	Φ 3.2 mm × 1. 1mm
9	Packing list	/
10	Instruction manual	/
11	Warranty card	/
12	Qualified certificate	/
13	Endochuck	/
14	Endo wrench	/
15	Water bottle	/
16	peristaltic pump	38 mm × 20 mm
17	Sealed ring	Φ 8.4 mm × 2 mm

The scaling tips and their accessories are not listed in this instruction manual completely. The detail can be found in the instruction for tips and packing list attached to the machine.

1.2.2 - Product performance and structure

Ultrasonic piezo scaler is composed of electrocircuit, waterway and ultrasonic transducer.

1.2.3 - Scope of application

Ultrasonic piezo scaler UDS-E is used for the dental calculus elimination and root canal treatment.

1.3 - THE MAIN TECHNICAL SPECIFICATIONS

- a) Adapter Input: 100V to 240V~ 50Hz/60Hz 1.2A (Max)
- b) Main unit input: 30VDC 1.3A
- c) Output primary tip Vibration excursion: $\leq 100\mu\text{m}$
- d) Output half-excursion force: $< 2\text{N}$
- e) Output tip vibration frequency: $28\text{kHz} \pm 3\text{kHz}$
- f) Output power: 3W to 20W
- g) Main unit fuse: 250VT 1.6AL
- h) Adapter fuse: 250VT 2.0AL
- i) Water pressure: 0.1bar to 5bar (0.01MPa to 0.5MPa)
- j) Weight of main unit: 1.35kg
- k) Weight of adapter: 0.3kg
- l) Operating mode: Continuous operation
- m) Type of protection against electric shock: Class II
- n) Degree of protection against electric shock: Type BF applied part
- o) Degree of protection against harmful ingress of water: Ordinary equipment (IPX0) Protection degree against water (used on the foot switch): IPX1
- p) Degree of safety of application in the presence of a Flammable Anaesthetic Mixture with air or with Oxygen or Nitrous Oxide: Equipment not suitable for being used in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide

1.4 - INSTALLATION OF THE MAIN COMPONENTS

1.4.1 - Sketch map for installation and connection.

a) The front and back sketch map of the main unit.


Fig. 1


Fig. 2

b) Sketch map for connection of foot switch, adapter and main unit.


Fig. 3

c) Water way connection.

* Sketch map of water way installation and connection


Fig. 4

* Sketch map water bottle installation


Fig. 5

d) Sketch map for connection of detachable handpiece.


Fig. 6

- e) Sketch map for how to install tip and endochuck with wrench.


Fig. 7

2 - PRODUCT FUNCTION AND USAGE

2.1 - SCALING FUNCTION

2.1.1 - Usage

- a) Open the packing box, make sure that all the parts and accessories are complete according to the packing list. Take the main unit out of the box and put it on a stable plane.
- b) Turn the water control switch to the max based on symbol as shown in 3.5.2 [note 1].
- c) Insert the plug of the foot switch to its socket(See Figure 3).
- d) Water way installation:
 - *Auto-water supply
Install the bottle reversely to the connector with adequate purified water (See Figure 5).
 - *Outside water
Connect one end of the water pipe to the water entrance, and the other end to the purified water source (See Figure 4).
- e) Screw the scaling tip tightly to handpiece by torque wrench, then connect the hand - piece and the connector of cable correctly.
- f) Insert the plug of the power source to its socket ,then get through to the power (See Figure 3).

- g) Switch on the main unit, then the scaling indicator and the fifth lead light of power regulator shines, auto-water indicator is on.
- h) According to need, push water button to choose auto-water or outside water. The indicator of choosen mode is on.
- i) Select a suitable scaling tip as you need, screw it on the handpiece tightly by the torque wrench (See Figure 7).
- j) The normal frequency is extremely high. Under the normal working state of scaling tips, a light touch and a certain to-and-fro motion will eliminate the tartar without heating. Overexertion and long-time lingering are forbidden.
- k) Vibrating intensity: Adjust the vibration intensity as you need, generally turn the knob to the middle grade. According to patient's different sensitivity and the rigidity of the gingival tartar, adjust the vibration intensity during the clinical treatment.
- l) Water volume adjust: Step on the foot switch, and the tip begins to vibrate, then turn the water control switch to form fine spray to cool down the handpiece and clean the teeth.
- m) The handpiece can be handled in the same gesture as a pen in hand.
- n) During the clinical treatment, be sure not to make the end of tip touch the teeth vertically and not to make the tip overexert on the surface of the teeth in case of hurting the teeth and damaging the tip.
- o) After finishing operation, keep the machine working for 30 seconds on the water supply condition in order to clean the handpiece and the scaling tip.
- p) Unscrew the scaling tip and pull out handpiece, then sterilize them.

Notice: Don't pull out the handpiece when the foot switch is stepped on and the machine is producing ultrasonic vibration.

2.1.2 - Instruction for main components of detachable handpiece (See Figure 6)

- a) Nipple: the nipple can be removed. You can screw out the nipple and clean the pole with alcohol termly.
- b) Handpiece seal: the seal can't be removed.
- c) Handpiece: The main part of the whole handpiece, can be autoclaved under the high temperature and pressure.
- d) The connector of the cable: connect the handpiece with the water source and power supply of the main unit.

Notice: Keep dry when the detachable handpiece connects to the connector of the cable.

2.1.3 Instruction of torque wrench (See Figure 7)

- a) The torque wrench's structure is designed in special way which can control the strength of the scaling tip's installation properly and correctly. It also can guarantee the operator screw or unscrew the scaling tip effectively and keep their hands away from being scratched.

b) Operation

- 1 Take the scaling tip into the torque wrench, operate as showed in figure 7.
- 2 Tip installation: Hold the handpiece, turn the tip toward direction as showed in figure 7 with the torque wrench. Turn two more circles when the tip stops, then the tip is installed.
- 3 Tip's uninstillation: Hold the handpiece, turn the wrench toward anti-clockwise direction.
- 4 Sterilize it in sterilizer after each treatment.
- 5 The torque wrench must be cooled naturally after sterilization to avoid scald when using next time.
- 6 Keep the torque wrench in a cool, dry and ventilated place and keep it clean.

2.2 - AUTO-WATER SUPPLY FUNCTION

a) Usage process

- 1 Pull out the water bottle vertically.
- 2 Open the lid, fill the bottle with adequate purified water and tighten the lid.
- 3 Clean the joint and interface of water bottle.
- 4 Place the bottle with water upside down, insert the joint into the interface vertically (See Figure 5).
- 5 Press "WATER" to choose auto-water supply.

b) **Notice**

- 1 Make sure the air hole and water entrance are not blocked. (See Figure 8).
- 2 Check if the cushion inside the lid is on good condition. If it becomes deformed or falls off, reinstall or change it.
- 3 Tighten the lid.
- 4 Clean the joint and interface of water bottle before each clinical operation.
- 5 After every clinical operation with clinical liquid, change a bottle with purified water, turn the water supply to max, make the machine work with auto-water supply for 30 seconds in order to keep the water way and spare parts clean and durable.
- 6 When the water is lower than water level lower limit, add water to keep water pipe smooth. (See Figure 9).


Fig. 8

Fig. 9

2.3 - ENDO FUNCTION

a) Usage process

- 1 Fix endochuck to handpiece by endo wrench (See Figure 7).
- 2 Unscrew the screw cap on the endochuck.
- 3 Put the ultrasonic file into the hole in the front of endochuck.
- 4 Screw down the screw cap with endo wrench to tight up the ultrasonic file.
- 5 Press option key, turn to endo function.
- 6 When ultrasonic scaler turns into endo function, only the first lead light is on and the power is at the 1st grade. Put the ultrasonic file into the patient's root canal slowly, step on the foot switch, During the treatment, turn up the power gradually according to the needs.

Notice

- 1 When fixing endochuck, it must be screwed down.
- 2 The screw cap on the endochuck must be screwed down.
- 3 Don't press it too much when the ultrasonic file in root canal.
- 4 Don't step on the foot switch until the ultrasonic file is in root canal.
- 5 The power range is supposed from the 1st to the 5th grades.

3 - STERILIZATION AND MAINTENANCE

3.1 - STERILIZATION OF DETACHABLE HANDPIECE

3.1.1 - Autoclaved in high temperature/pressure:

- a) 121°C/1bar (0.1MPa), 20 minutes.
- b) 135°C/2.2bar (0.22MPa), 18 minutes.
- c) Pull out the handpiece and unscrew scaling tip and endochuck after each operatin.
- d) Pack the handpiece with sterile gauze or sterile bag before sterilizing.
- e) Reuse handpiece after it cools naturally in case of scaling hand.

3.1.2 - Notice

- a) Dry the cleaning liquid in the handpiece with compressed air before sterilization.
- b) Be sure that the scaling tip has been unscrewed from the handpiece and it cannot sterilized with others.
- c) Please notice whether the outer of the handpiece is damaged during the treatment or sterilization, don't smear any protective oil on the surface of handpiece.
- d) There are two waterproof "o" rings at the end of handpiece. Please lubricate them with dental lube frequently, as sterilization and repeated pulling and inserting will reduce their life. Change a new one once it is damaged or worn excessively.
- e) The following sterilizing methods are forbidden:
 - 1 Put handpiece into any liquid for boiling.
 - 2 Dip handpiece in disinfectors such as iodine, alcohol and glutaraldehyde.
 - 3 Put handpiece into oven or microwave oven for baking.

3.2 - STERILIZATION OF SCALING TIPS AND ENDOCHUCK

- a) All the scaling tips and endochuck can be autoclaved under high temperature and pressure.

3.3 - STERILIZATION OF TORQUE WRENCH AND ENDO WRENCH

- a) The torque wrench and endo wrench can be sterilized in high temperature and pressure.
- b) The following sterilization ways for torque wrench are forbidden:
- 1 Braise in liquor.
 - 2 Dip in iodine, alcohol or glutaraldehyde.
 - 3 Torrefy in oven or microwave oven.

Notice: We are not responsible for any damage of the torque wrench directly or indirectly made by any way in the above items.

3.4 - CLEANING OF TIPS, ENDOCHUCK, TORQUE WRENCH AND ENDO WRENCH

The scaling tip, endochuck, torque wrench and endo wrench can be cleaned by ultrasonic cleaner.

3.5 - TROUBLESHOOTING AND NOTES

3.5.1 - Troubleshooting

Fault	Possible cause	Solutions
The scaling tip doesn't vibrate and there is no water flowing out when stepping on the foot switch.	The power pipe plug is in loose contact.	Make the plug insert to the socket well.
	The foot switch is in loose contact.	Insert the foot switch to its socket tightly.
	The fuse in the main unit is broken.	Contact our dealers or us.
The scaling tip doesn't vibrate but there is water flowing out when stepping on the foot switch.	The tip is in loose contact.	Screw the tip on the handpiece tightly (See Figure 7).
	The connect plug between the handpiece and the circuit board is in loose contact.	Contact our dealers or us.
	Something wrong with the handpiece.	Send the handpiece to our company to repair.
	Something wrong with the cable.	Contact our dealers or us.
The scaling tip vibrates but there is no spray when stepping on the foot switch.	The water control switch is not on.	Turn on the water control switch [note 1].
	There is impurity in the electric-magnetic valve.	Contact our dealers or us.
	The water system is blocked.	Clean the water pipe by multi-function syringe [note 2].
There is still water flowing out after the power is off.	There is impurity in the electric-magnetic valve.	Contact our dealers or us.
The handpiece generates heat.	The water control switch is in a low setting.	Turn the water control switch to a higher grade [note 2].

Fault	Possible cause	Solutions
The amount of spouting water is too little.	The water pressure is not high enough.	Make the water pressure higher.
	The water pipe is blocked.	Clean the water pipe by multi-function syringe [note 2].
The vibration of the tip becomes weak.	The tip hasn't been screwed on to the handpiece tightly.	Screw the tip on the handpiece tightly (See Figure 7).
	The tip is loose because of vibration.	Screw on the tip tightly (See Figure 7).
	The coupling between the handpiece and the cable isn't d dry.	Dry it by the hot air.
	The tip is damaged [note 3].	Change a new one.
There is water seeping from the coupling between the handpiece and the cable.	The waterproof "O " ring is damaged.	Change a new "O " ring.
The u-file doesn't vibrate.	The screw is loose.	Tighten it.
	Endochuck is damaged.	Change a new one.
There is noise coming from the endochuck.	The screw is loose.	Tighten it.
Peristaltic Leaks.	The inner water pipe cracks.	Replace a new peristaltic pump.
There is no water coming out from the handpiece (automatic water supply mode).	There is air in the water pipe.	Turn the water control to the Max, reinsert the bottle.

If the problem still can't be solved, please contact with local dealer or manufacturer.

3.5.2 - Notes

- a) [Note 1] Turn the water control switch toward direction showed in symbol till not be able to turn anymore, it comes to the min, on the contrary direction, the water volume increases.
- b) [Note 2] Clean the water pipe with the Multi-function syringe of the dental unit (See Figure 10):


Fig. 10

- 1 Cut the water pipe at a distance of 10cm to 20cm from the water entrance.
 - 2 Turn on the electricity and get through to the power.
 - 3 Connect the multi-function syringe of dental unit to the water pipe.
 - 4 Disassemble the tip or handpiece.
 - 5 Step on the foot switch.
 - 6 Turn on the switch of the multi-function syringe, press the water into the machine and the impurity blocked in the water pipe can be eliminated.
- c) [Note 3] If the scaling tip has been screwed on tightly and there is fine spray too, the following phenomena show that the scaling tip is damaged:
- 1 The vibrating intensity and the water spray degree become weak obviously.
 - 2 During treatment, it produces the sound like "buzz" from the scaling tip.

d) Note 4.

- 1 Pull out the water bottle vertically.
- 2 Remove the bolt with screwdriver. (See Figure 11)
- 3 Pull the two fasten rings out to the right, then pull the peristaltic pipe out to the right. (See Figure 12)
- 4 Pull the peristaltic out from the drive pole, and install the new peristaltic back to the same position. (See Figure 12)
- 5 Install the fasten rings onto the ends of pipes then insert the pipes into water-in and water-out joints. Push the fasten rings to mid of water-in and water-out joints. (See Figure 12)
- 6 Install the lid to the right position.


Fig. 11

**Notice:**

- 1 Dry the drive pole to avoid being slippery.
- 2 Install the peristaltic onto the drive pole carefully and lightly to avoid the damage of the inner structure of the machine.
- 3 After changing a new peristaltic, turn water supply to max and make machine work with auto-water supply for 30 seconds in order to be compatible with peristaltic and its pipes.


Fig. 12

4 - PRECAUTION

4.1 - NOTICE WHEN USING EQUIPMENT

- Keep the scaler clean before and after operation.
- The handpiece, scaling tip, torque wrench, endo wrench and endochuck must be sterilized before each treatment.
- Don't screw or unscrew the scaling tip and endochuck when stepping on the foot switch.
- The scaling tip must be fastened and there must be fine spray or drip coming out from the tip when operating.
- Change a new one when the tip and ultrasonic file are damaged or worn excessively.
- Don't twist the tip and endochuck or rub them.
- If use the water source without hydraulic pressure, the water surface should be one meter higher than the head of the patient.
- Insure the connector of handpiece and the socket of the cable dry before installing the handpiece.
- Don't pull the cable forcibly in case of the handpiece falling off the cable.
- Don't knock or rub the handpiece.
- After operating, turn off the power, then pull out the plug.
- We are only responsible for the safety on the following conditions:
 - I The maintenance, repairment and modification are made by the manufacturer or the authorized dealer.
 - II The changed components are original of "Carlo De Giorgi" and are operated according to instruction manual.
- The internal screw thread of the scaling tips produced by some manufacturers is coarse, rusty and collapsed. This will damage the external screw thread of the handpiece irretrievably. Please use "Carlo De Giorgi" brand scaling tips.
- This model is only matched the adapter of our company.
- Keep peristaltic dry. Otherwise, its wheels are likely to be slippery and the water supply will not work.

4.2 - CONTRAINDICATION

- The hemophilia disease patient is not allowed to use this equipment.
- The patients or doctors with heart pacemaker are forbidden to use this equipment.
- The heart disease patient, pregnant woman and children should be cautious to use the equipment.

4.3 - STORAGE AND MAINTENANCE

- The equipment should be handled carefully and lightly. Be sure that it is far from the vibration, and installed or kept in a cool, dry and ventilated place.
- Don't store the machine together with the articles that are combustible, poisonous, caustic, and explosive.
- This equipment should be stored in a room where the relative humidity is $\leq 80\%$, atmospheric pressure is 50kPa to 106kPa, and the temperature is -10°C to $+50^{\circ}\text{C}$.
- If not use for a long time, please make the machine get through to the power and water once per month for five minutes.

4.4 - TRANSPORTATION

- Excessive impact and shake should be prevented in transportation. Lay it carefully and lightly and don't invert it.
- Don't put it together with dangerous goods during transportation.
- Avoid solarization and getting wet in rain and snow during transportation.

4.5 - WORKING CONDITION

- a) Environment temperature: 5°C to 40°C .
- b) Relative humidity: $\leq 80\%$.
- c) Atmosphere pressure: 50kPa to 106kPa.



5 - AFTER SERVICE

We offer one year warranty.

The repair of the equipment should be carried out by our professional technician. We are not responsible for any irretrievable damage caused by the non-professional person.

6 - SYMBOLS INSTRUCTION

	Trademark.		Caution, consult accompanying documents.
	Class II equipment.		Use indoor only.
	Type BF applied part.		Can be autoclaved.
	Alternating current.		30 VDC power supply socket.
	Connection for the foot switch.		Water entrance pressure 0.01 Mpa-0.5MPa.
	Adjustment for the water flow.		FDA marked product.
	Temperature limitation.		CE marked product.
	Humidity limitation.		Atmospheric pressure limitation.
	Date of manufacture.		Appliance compliance WEEE directive.
	Authorized Representative in the EUROPEAN COMMUNITY.		Manufacturer.
	Add water.		Mode of outside water source.
	Mode of auto-water system.		

7 - ENVIRONMENTAL PROTECTION

There is not any harmful factor in our product. You can deal with it based on the local law.

We reserve the rights to change the design of the equipment, the technique, fittings, the instruction manual and the content of the original packing list at any time without notice. If there are some differences between blueprint and real equipment, take the real equipment as the norm.

8 - DECLARATION OF CONFORMITY

8.1 - PRODUCT CONFORMITY THE FOLLOWING STANDARDS:

EN ISO 14971:2000+A1:2003
EN 980: 2003 / ISO 9687:1993 / EN 1041:1998
EN 60601-1:1990+A1:1993+A2:1995+ A13:1996
EN 60601-1-4:2001 / EN 60601-1-2:2001 + A1:2006
EN ISO 7405:1997
EN ISO 10993-1:2003 / EN ISO 10993-5:1999
EN ISO 10993-10:2003+A1:2003
EN 61205:1994 / EN ISO 22374:2005
EN ISO 17664:2004 / EN ISO 17665-1:2006

9 - FOR TECHNICAL DATA, PLEASE CONTACT

Carlo De Giorgi Srl
Via Tonale 1/3
20021 Baranzate
Milano Italy
Tel. 02 356 15 43
Fax 02 356 18 08
<http://www.degiorgi.it>

9.1 - EMC - DECLARATION OF CONFORMITY

Guidance and manufacturer's declaration - electromagnetic emission		
The model UDS-E is intended for use in the electromagnetic environment specified below. The customer or the user of the model UDS-E should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The model UDS-E uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR11	Class B	The model UDS-E is suitable for used in domestic establishment and in establishment directly connected to a low voltage power supply network which supplies buildings used for domestic purposes.
Harmonic emission IEC 61000-3-2	ns Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	

**Guidance and manufacturer's declaration - electromagnetic emission**

The model UDS-E is intended for use in the electromagnetic environment specified below. The customer or the user of the model UDS-E should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment -
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 % .
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1 kV for Input/output lines	±2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11.	< 5 % UT (> 95 % dip in UT.) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70% UT (30 % dip in UT) for 25 cycles < 5 % UT (> 95 % dip in UT) for 5 sec	< 5 % UT (> 95 % dip in UT.) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70% UT (30 % dip in UT) for 25 cycles < 5 % UT (> 95 % dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment . If the user of the model UDS-E requires continued operation during power mains interruptions , it is recommended that the model UDS-E should be powered from an uninterruptible power supply or a battery.
Power frequency (50 / 60 Hz) magnetic field IEC 61000-4-8	3 A/m	Not applicable	Not applicable

NOTE UT is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic emission			
The model UDS-E is intended for use in the electromagnetic environment specified below. The customer or the user of the model UDS-E should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	150 kHz to 80 MHz	3V	Portable and mobile RF communications equipment should be used no closer to any part of the model UDS-E, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $3V$
Radiated RF IEC 61000-4-3	3 Vrms 80 MHz to 2.5 Ghz	3 V/m	3 V/m $d=1.2 \times P^{1/2}$ 80 MHz to 800 MHz $d=2.3 \times P^{1/2}$ 800 MHz to 2.5 Ghz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz end 800 MHz. the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model UDS-E is used exceeds the applicable RF compliance level above, the model UDS-E should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the model UDS-E.

**Recommended separation distances between
portable and mobile RF communications equipment and the model
UDS-E**

The model UDS-E is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the model UDS-E can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the model UDS-E as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80MHz $d = 1.2 \times P^{1/2}$	80MHz to 800MHz $d = 1.2 \times P^{1/2}$	800MHz to 2,5GHz $d = 2.3 \times P^{1/2}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) accordable to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz. the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

The device has been tested and homologated in accordance with EN 60601-1-2 for EMC. This does not guarantee in any way that this device will not be effected by electromagnetic interference. Avoid using the device in high electromagnetic environment.

carlo de giorgi s.r.l.

Via Tonale, 1 – 20021 Baranzate (MI) – Italy
Tel. +39.02.356.15.43 r.a. – Fax +39.02.356.18.08
<http://www.degiorgi.it> E-Mail: info@degiorgi.it

MAN 032 REV 01/2010