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	UDS-E (LED) UDS-L (LED) UDS-A (LED) UDS-P (LED)	UDS-J2 (LED)	UDS-J UDS-B UDS-K (LED)	UDS-N3 (LED)	UDS-N1 UDS-N2 (LED) UDS-N4	Marca di ablatore compatibile
Punte da Ablatore						
G1	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	Compatibile con ablatori Woodpecker ed EMS
G2	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	
G3	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	
G4	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	
G5	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	
G6	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	
G7	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	
G8	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	
G9	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	
G10	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	
G11	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	
Punte da Perio						
P1	1-10 (P)	1-7	1-6	LOW-MID	LOW-MID	Compatibile con ablatori Woodpecker ed EMS
P2L	1-3(P)	1-2	1-2	LOW	LOW	
P2LD	1-2(P)	1	1	LOW	LOW	
P2R	1-3(P)	1-2	1-2	LOW	LOW	
P2RD	1-2(P)	1	1	LOW	LOW	
P3	1-6(P)	1-4	1-3	LOW-MID	LOW-MID	
P3D	1-6(P)	1-4	1-3	LOW-MID	LOW-MID	
P4	1-6(P)	1-4	1-3	LOW-MID	LOW-MID	
Punte da Endo						
E1	1-3(E)	-	-	LOW	-	Compatibile con ablatori Woodpecker ed EMS
E2	1-3(E)	-	-	LOW	-	
E3	1-6(E)	-	-	LOW	-	
E3D	1-3(E)	-	-	LOW	-	
E4	1-6(E)	-	-	LOW	-	
E4D	1-3(E)	-	-	LOW	-	
E5	1-6(E)	-	-	LOW	-	
E5D	1-3(E)	-	-	LOW	-	
P4D	1-6(E)	-	-	LOW	-	
E8	1-10(E)	-	-	LOW	-	
E9	1-10(E)	-	-	LOW	-	
E10	1-6(E)	-	-	LOW	-	
E10D	1-6(E)	-	-	LOW	-	
E11	1-6(E)	-	-	LOW	-	
E11D	1-6(E)	-	-	LOW	-	
E14	1-3(E)	-	-	LOW	-	
E15	1-3(E)	-	-	LOW	-	
Preparazione delle carie						
SB1	1-10 (P)	1-7	1-6	LOW-MID	LOW-MID	Compatibile con ablatori Woodpecker ed EMS
SB2	1-10 (P)	1-7	1-6	LOW-MID	LOW-MID	
SB3	1-10 (P)	1-7	1-6	LOW-MID	LOW-MID	
SBL	1-10 (P)	1-7	1-6	LOW-MID	LOW-MID	
SBR	1-10 (P)	1-7	1-6	LOW-MID	LOW-MID	

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Tel. +39.02.356.15.43 r.a. – Fax +39.02.356.18.08

CONTENTS

1 - DEVICE COMPONENTS AND SETUP	4
1.1 - INSTRUCTIONS.....	4
1.2 - COMPONENTS.....	5
1.2.1 – COMPONENTS INCLUDED	5
1.2.2 - PRODUCT PERFORMANCE AND STRUCTURE.....	5
1.2.3 - SCOPE OF APPLICATION	5
1.3 - MAIN TECHNICAL SPECIFICATIONS	6
1.3.1 - TECHNICAL SPECIFICATIONS OF THE ULTRASONIC SCALER	6
1.4 - INSTALLATION OF THE MAIN COMPONENTS	7
1.4.1 - PRODUCT COMPONENTS	7
1.4.2 – HANDPIECE COMPONENTS.....	8
1.4.3 - TORQUE WRENCH USE INSTRUCTIONS	9
2. - PRODUCT FUNCTION AND USE	11
2.1 - OPERATION	11
2.2 - OPERATION METHOD AND FUNCTION INSTRUCTION	12
3 - STERILIZATION AND MAINTENANCE	13
3.1 –STERILIZATION OF THE HANDPIECE	13
3.1.1 - ACCEPTABLE AUTOCLAVE TEMPERATURE/PRESSURE.....	13
3.1.2 - NOTICE	13
3.2 - STERILIZATION OF THE SCALING TIP.....	14
3.3 - STERILIZATION OF THE TORQUE WRENCH	14
3.4 - CLEANING OF THE TORQUE WRENCH AND SCALING TIPS..	14
3.5 - STERILIZATION AND CLEANING OF THE LED LAMP AND LAMP PIPE	14
3.5 – TROUBLESHOOTING AND NOTES.....	15
4 - PRECAUTIONS.....	18
4.1 - NOTICE	18
4.2 - CONTRAINDICATIONS	20
4.3 - STORAGE AND MAINTENANCE	20
4.4 - TRANSPORTATION	20
5 - AFTER-SERVICE	21
6 - SYMBOL INSTRUCTION.....	21
7 - ENVIRONMENTAL PROTECTION.....	23

8 - MANUFACTURER’S RIGHT 23

9 - AUTHORIZED REPRESENTATIVE IN THE EU..... 23

10 – DECLARATION OF CONFORMITY 24

11 – STATEMENT 29

TABLE OF OPERATING POWER OF THE TIPS..... 30

1 - DEVICE COMPONENTS AND SETUP

1.1 - INSTRUCTIONS

The UDS-K ultrasonic piezo scaler is designed primarily for teeth cleaning, and it's also an indispensable equipment for tooth disease prevention and treatment. It contains the following features:

1. Automatic frequency tracking, to ensure that the machine always works steadily on the best frequency.
2. The handpiece is detachable and can be autoclaved at a high temperature of 135°C and a pressure of 0.22MPa.
3. Digital controlled, easy operation and more efficient for scaling.

1.2 - COMPONENTS

1.2.1 – Components included

The following parts are included in a complete UDS-K/UDS-K LED ultrasonic piezo scaler:

No.	UDS-K	UDS-K LED
01	Main unit	Main unit
02	Detachable handpiece	Detachable handpiece
03	Power supply	Power supply
04	Foot pedal	Foot pedal
05	Scaling tips	Scaling tips
06	Water line	Water line
07	Torque wrench	Torque wrench
08	Waterproof “O” ring	Waterproof “O” ring
09	Handpiece seal	Handpiece seal
10	Packing list	Packing list
11	Instruction manual	Instruction manual
12	-	LED light

1.2.2 - Product performance and structure

The ultrasonic scaler is composed of electrocircuit, waterway and ultrasonic transducer.

1.2.3 - Scope of application

The ultrasonic piezo scaler UDS-K is used for dental calculus elimination.


1.4 - INSTALLATION OF THE MAIN COMPONENTS


1.4.1 - Product components

The components of the equipment are shown in picture 1:

1.4.2 – Handpiece components

The components of the detachable handpiece are shown in picture 2:

Picture 2

Instruction for main components of detachable handpiece.

Nipple: the nipple can be removed. You can screw out the nipple and clean the shaft with alcohol.

Decorative ring: the ring can be removed and cleaned with alcohol regularly.

Handpiece: the main body can be autoclaved under high temperature and pressure.

Symbol: autoclavable (135°C, 0.22 MPa).


Cable connector: connects the handpiece with the water source and power supply of the main unit.

Note: keep the connection points dry.

1.4.3 - Torque wrench use instructions

a) Brief introduction and description

The structure of the torque wrench is designed specifically to control the strength of the scaling tip installation properly and correctly. It also avoids damage to the operator hands.


b) Sterilization

Sterilize in steam with a temperature of 135°C, and a pressure of 0.22 MPa).

c) Operation

1. Use the wrench as shown in picture 4.

Picture 4

2. Tip installation. Hold the handpiece, turn the tip in the clockwise direction with the torque wrench. Turn two more circles after the tip stops, then the tip is installed.
3. Tip removal. Hold the handpiece, turn torque wrench in the counter-clockwise direction.
4. Sterilize after each treatment.
5. The torque wrench must be cooled naturally after sterilization before using again, to avoid scalding.
6. Keep the torque wrench clean, in a cool, dry and ventilated place.

d) Warning.

The following sterilization methods are forbidden:

1. Dipping in liquor;
2. Dipping in iodine, alcohol or glutaraldehyde;
3. Heat sterilization in oven or microwave oven.

Note: we are not responsible for the damage of the torque wrench for any case listed in above.

2. - PRODUCT FUNCTION AND USE

2.1 - OPERATION

- a) Open the packing box, make sure that all the parts and accessories are complete according to the components list.
- b) Take the main unit of the box and put it on a stable flat surface.
- c) Turn the water control switch all the way clockwise and turn the vibrating intensity control knob clockwise in a suitable position.
- d) Insert the plug of the foot pedal in its socket.
- e) Connect one end of the water pipe to the water entrance, and the other end to a clean water source.
- f) Connect the handpiece: screw on the scaling tip to the handpiece with the torque wrench, then insert the connector of the cable to the handpiece correctly.
- g) Connect the device to the power outlet.
- h) Press the power switch of the main unit, and verify that the power indicator turns on.

2.2 - OPERATION METHOD AND FUNCTION INSTRUCTION

- a) Make the scaler easily accessible to the operator. Before turning on, please turn the vibrating intensity control knob to the minimum, and the water control switch to the maximum (turn three circles towards clockwise direction from the minimum to the maximum).
- b) The normal frequency is as high as 28KHz $\pm$ 3kHz. Under normal working conditions, a light touch and a certain to-and-fro motion will eliminate the tartar without excessive heating. Overexertion and overstay are forbidden.
- c) Select a suitable scaling tip according to the treatment, screw on it to the handpiece tightly with the torque wrench (as shown in picture 4).
- d) Step on the foot switch, the tip begins to vibrate. In the UDS-K LED, the LED light on the top of the handpiece also turns on. After releasing the foot switch, the lights will stay on for 10 more seconds.
- e) Vibrating intensity: Adjust the vibrating intensity as needed; generally turn the knob to the middle grade. Because different patients have different sensitivity and different rigidity of the gingival tartar, the vibrating intensity should be adjusted during the clinical treatment.
- f) Water adjust: step on the foot pedal, the tip begins to vibrate, then turn the water control switch to form fine spray to cool down the handpiece and clean the teeth.
- g) The handpiece can be handled in the same gesture as a pen in hand.
- h) Let the tip touch the surface of the teeth lightly; do not use too much pressure, or else the teeth will be hurt and the scaling tip will be damaged.
- i) After the treatment, keep the machine working for 30 seconds with the water supply open, so that the handpiece and the scaling tips can be cleaned.
- j) Remove the handpiece and unscrew the scaling tip, and sterilize both.

Please note: Don't pull out the handpiece when the foot pedal is pressed and the machine is working.

3 - STERILIZATION AND MAINTENANCE

3.1 –STERILIZATION OF THE HANDPIECE

3.1.1 - Acceptable autoclave temperature/pressure

- a) Sterilization at 121°C/1 Bar (0,1 MPa);
- b) Sterilization at 135°C/2.2 Bar (0,22 MPa).

Remove the handpiece and unscrew the scaling tip after treatment.
Pack the handpiece with sterile gauze or sterilization bag.
Wait for the handpiece to cool naturally before use, to avoid scalding.

3.1.2 - Notice

- a) Clean the cleaning liquid in the handpiece with compressed air before sterilization.
- b) Be sure that the scaling tip has been unscrewed from the handpiece, before sterilization.
- c) Please pay attention whether the shell of the handpiece have been damaged during the treatment and sterilization. Don't smear any protective oil on the surface of handpiece.
- d) There are two waterproof "O" rings at the end of the handpiece. Please lubricate them with dental lube frequently, as sterilization and repeated pulling and inserting will reduce their use life. Once they are damaged or excessively worn, please change with new ones.
- e) The following sterilizing methods are forbidden:
 - 1. Put the handpiece into any liquid for boiling.
 - 2. Put the handpiece in disinfectors such as alcohol, iodine or glutaraldehyde.
 - 3. Put the handpiece in oven or microwave oven for baking.

3.2 - STERILIZATION OF THE SCALING TIP

All the scaling tips can be autoclaved at 135°C.

3.3 - STERILIZATION OF THE TORQUE WRENCH

The wrench can be sterilized under high temperature and pressure.

3.4 - CLEANING OF THE TORQUE WRENCH AND SCALING TIPS

Both of them can be cleaned with an ultrasonic cleaner.

3.5 - STERILIZATION AND CLEANING OF THE LED LAMP AND LAMP PIPE

Clean the LED lamp and lamp pipe with purified water and sterilize them under high temperature and pressure after every treatment.

3.5 – TROUBLESHOOTING AND NOTES

Fault	Possible causes	Solutions
The scaling tip does not vibrate and there is no water flowing out when getting through the electricity.	The power plug is in loose contact.	Make the plug insert to the socket well.
	The foot switch is in loose contact.	Insert the foot switch in its socket tightly.
	The fuse of the transformer is broken.	Contact our dealers or us.
	The fuse in the main unit is broken.	Contact our dealers or us.
The scaling tip does not vibrate but there is water flowing out when pressing the foot switch.	The handpiece is damaged.	Contact our dealers or us.
	The cable is damaged.	Contact our dealers or us.
The scaling tip vibrates but there is no water flowing out when pressing the foot switch.	The water control switch is turned off.	Turn on the water control switch [note 1].
	There is impurity in the solenoid valve.	Contact our dealers or us.
	The solenoid valve is damaged.	Contact our dealers or us.
There is still water flowing out after the power is off.	There is impurity in the electric-magnetic valve.	Contact our dealers or us.
The handpiece generates heat.	The water control switch is in a low grade.	Turn the water control switch to a higher grade [note 1].
The amount of water output is too little.	The water control switch is in a low grade.	Turn the water control switch to a higher grade [note 1].
	The water pressure is not high enough.	Turn the water control switch to a higher grade [note 1].
	The water pipe is blocked.	Clean the water line by multi-function syringe [note 2].


Fault	Possible causes	Solutions
The vibration of the tip becomes weak.	The tip hasn't been screwed on to the handpiece tightly.	Screw on the tip to the handpiece tightly (as showed in the picture 4).
	The tip was shaken loose by vibration	Screw on the tip tightly (as showed in the picture 4).
	The coupling between the handpiece and the cable isn't dry.	Dry it with hot air.
	The tip is damaged [note 3].	Change with a new one.
The vibration is too strong and the vibrating intensity control knob is malfunctioning.	The vibration intensity control knob is damaged.	Contact our dealers or us.
There is water seeping from the coupling between the handpiece and the cable.	The waterproof "O" ring is damaged.	Change with a new "O" ring.

If the troubles still can't be solved yet, please contact the local distributors or our company.

[Note 1] The water control knob can adjust the water volume according to the symbol.

[Note 2] Clean the water pipe with the multi-function syringe of the dental unit (as showed in the picture 5):

- Snip the water pipe with scissors at a distance of 10cm-15cm from the water entrance.
- Turn on the power switch.
- Connect the multi-function syringe of the dental unit to the water pipe.
- Press the foot switch.
- Turn on the switch of the three-way syringe, press the air or water into the water pipe in the machine, then the impurity blocked in the water line can be eliminated.

[Note 3] If the scaling tip has been screwed on tightly and the water output is correct, the following phenomena show that the scaling tip is damaged:

- The vibrating intensity and the water atomization decrease evidently.
- During treatment, there is a significant buzzing sound coming from the scaling tip.

4 - PRECAUTIONS

4.1 - NOTICE

- Keep the scaler clean before and after operation.
- The handpiece, scaling tip, and torque wrench must be sterilized before every treatment.
- Don't screw or unscrew the scaling tip while pressing down on the foot pedal.
- The scaling tip must be tightly screwed and there must be fine spray coming out from the tip when operating.
- Replace with a new one when the tip is damaged or excessively worn.
- During the treatment, the tip will heat up considerably if the water flow is too low. Make sure the water flow is adequate.
- Don't twist or file the tip.
- Don't use an impure water source. Never replace the distilled water with physiological solution.
- If using a water source without hydraulic pressure, the water surface should be one meter higher than the head of the patient.
- Make sure that the connecting end of handpiece and the socket of the cable connector are completely dry before handpiece installation.
- Don't pull the cable hard during treatment.
- Avoid shocks or hits to the handpiece.
- After treatment, turn off the power, then pull out the plug.
- Plug the device in an easily accessible power socket, to make sure it can be unplugged easily in case of an emergency.
- The power supply is considered integral part of the device. Do not use the device with a different power supply.
- The power supply is not waterproof. Keep it dry and away from water.
- As a professional manufacturer producing medical instruments, we are responsible for the safety only when maintenance, repairs and any changes are done by Carlo De Giorgi s.r.l. or our authorized distributors. Only original accessories should be used, and the device should be operated by the manual.

-
- The screw thread scaling tips produced by some other manufacturers may be coarse, rusty and damaged; this will damage the screw thread of the handpiece irretrievably. Please use only original scaling tips.
 - Select a suitable power level when using different types of tips (see "TABLE OF OPERATING POWER OF THE TIPS").

4.2 - CONTRAINDICATIONS

- Do not use on patients with haemophilia.
- Not for use by doctors, or on patients, with an artificial pacemaker.
- Pay extra care when using the device on patients with heart disease, pregnant women, and children.

4.3 - STORAGE AND MAINTENANCE

- The device should be handled carefully and lightly. Be sure that it is far from external vibrations, and install or keep it in a cool, dry and ventilated place.
- Don't store the device with combustible, poisonous, caustic, or explosive material.
- This device should be stored in a room where the relative humidity is $\leq 80\%$, atmospheric pressure is 50kPa ~ 106kPa, and the temperature is $-10\text{ }^{\circ}\text{C} \sim +50\text{ }^{\circ}\text{C}$.
- Please turn off the power switch and pull out the power line plug when the device is not used. If the device is used for a long time, please connect it to the power line and let the water flow through for five minutes once a month.

4.4 - TRANSPORTATION

- Excessive impacts and hits should be avoided during transportation. Lay the device carefully and lightly and don't overturn it.
- Don't store the device together with dangerous goods.
- Avoid solarisation and getting wet in rain and snow during transportation.

5 - AFTER-SERVICE

- We offer free repairs for any production defect of the device for one year's from the date when it is sold to the final customer. We are not responsible for damage to the device deriving from improper use.
- Any repair of the device should be carried out by our professional technicians. We are not responsible for any damage caused by external technicians.

6 - SYMBOL INSTRUCTION

Manufacturer trademark

CE marked product

Date of manufacture

Manufacturer

Insulation type BF

Insulation class II

IPX1

Anti-drip device

Recycling

Use indoors

Keep dry

Autoclavable

Fragile, handle with care


Alternating current

Foot switch interface


Water flow setting


Storage atmospheric pressure

H₂O
0.01-0.5MPa

Water inlet pressure: 0.01MPa - 0.5MPa

Acceptable moisture limits

Acceptable temperature limits

WEEE Directive implementation

Please read the user manual

Authorised Representative in the EUROPEAN COMMUNITY

7 - ENVIRONMENTAL PROTECTION

There are no harmful factors in our product. You can dispose with it based on the local law.

8 - MANUFACTURER'S RIGHT

We reserve the right to change the design of the equipment, the technique, fittings, the instruction manual and the content of the original packing list at any time without notice. If there are some differences between blueprint and real equipment, take the real equipment as the norm.

9 - AUTHORIZED REPRESENTATIVE IN THE EU

Wellkang Ltd (www.CE-Marking.eu)
29 Harley St., LONDON, W1G 9QR, UK

10 – DECLARATION OF CONFORMITY

桂林市啄木鸟医疗器械有限公司
GUILIN WOODPECKER Medical Instrument Co.,LTD.

File No.: ZMN/WI-04-031-02

Version: F

EC Declaration of Conformity

Manufacturer:

whose single Authorized Representative:

Guilin Woodpecker Medical Instrument Co., Ltd.

Wellkang Ltd t/a Wellkang Tech Consulting

We, the manufacturer, herewith declare that
the products **ULTRASONIC SCALER**, UMDNS-Code: 17497

MODEL: UDS-K

meet the provisions of Directive 93/42/EEC which apply to them.

The medical device has been assigned to class IIa according to Annex IX of the Directive 93/42/EEC. It bears the mark

CE 0197

The product concerned has been manufactured under a quality management system according to Annex V of Directive 93/42/EEC.

Compliance of the designated product with the Directive 93/42/EEC has been assessed and certified by the Notified Body

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Issue date: 2015-08-17

Expiry date: 2019-09-15

following the procedure relating to the EC Declaration of Conformity set out in Annex V of Directive 93/42/EEC.

This Declaration of Conformity covers all medical devices as specified in the product list belonging to this declaration and is only valid in connection with a batch specific Certificate of Compliance for all products concerned bearing the CE mark

The above mentioned declaration of conformity is exclusively under the responsibility of

Company: Guilin Woodpecker Medical Instrument Co., Ltd.

Address: Information Industrial Park, Guilin National High-Tech Zone, Guilin, GuangXi,
541004, P.R.China

2015.8.17
Place, date

Legally binding signature, Function

桂林市啄木鸟医疗器械有限公司
GUILIN WOODPECKER Medical Instrument Co.,LTD.

File No.: ZMN/WI-04-131-02

Version: E

EC Declaration of Conformity

Manufacturer:

whose single Authorized Representative:

Guilin Woodpecker Medical Instrument Co., Ltd.

Wellkang Ltd t/a Wellkang Tech Consulting

We, the manufacturer, herewith declare that
the products **ULTRASONIC SCALER**, UMDNS-Code: **17497**

MODEL:UDS-K LED

meet the provisions of Directive 93/42/EEC which apply to them.

The medical device has been assigned to class IIa according to Annex IX of the Directive 93/42/EEC. It bears the mark

CE 0197

The product concerned has been manufactured under a quality management system according to Annex V of Directive 93/42/EEC.

Compliance of the designated product with the Directive 93/42/EEC has been assessed and certified by the Notified Body

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Issue date: 2015-08-17

Expiry date: 2019-09-15

following the procedure relating to the EC Declaration of Conformity set out in Annex V of Directive 93/42/EEC.

This Declaration of Conformity covers all medical devices as specified in the product list belonging to this declaration and is only valid in connection with a batch specific Certificate of Compliance for all products concerned bearing the CE mark

The above mentioned declaration of conformity is exclusively under the responsibility of

Company: Guilin Woodpecker Medical Instrument Co., Ltd.

Address: Information Industrial Park, Guilin National High-Tech Zone, Guilin, GuangXi,
541004, P.R.China

Place, date

2015.8.17

[Handwritten signature]

Legally binding signature, Function

**Guide & Manufacturer declaration - Electromagnetic emissions**

The device UDS-K\UDS-K LED is intended for use in the electromagnetic environment specified below. The customer or the user of the device UDS-K\UDS-K LED should make sure that it is used in such an

Emission test	Compliance	Electromagnetic environment - guide
RF emissions CISPR 11	Group 1	The device UDS-K\UDS-K LED uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The device UDS-K\UDS-K LED is suitable for used in domestic environments and in establishments directly connected to a low voltage power supply network which supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Guide & Manufacturer declaration - Electromagnetic immunity

The device UDS-K\UDS-K LED is intended for use in the electromagnetic environment specified below. The customer or the user of the device UDS-K\UDS-K LED should make sure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guide
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1 kV for input/output lines	±2kV for power supply lines ±1kV for interconnecting cable	The main power should be of the kind of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line to line ±2 kV line to ground	±1 kV line to line	The main power should be of the kind of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11.	<5 % U_T (>95% dip in U_T) for 0.5 cycle 40 % U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycle <5% U_T (>95 % dip in U_T) for 5 sec	<5 % U_T (>95% dip in U_T) for 0.5 cycle 40 % U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycle <5% U_T (>95 % dip in U_T) for 5 sec	The main power should be of the kind of a typical commercial or hospital environment. If the user of the device UDS-K\UDS-K LED requires a continuous operation during power supply interruptions, we recommend that the UDS-K\UDS-K LED be connected to an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) Magnetic field IEC 61000-4-8	3 A/m	3 A/m	The magnetic field should be of the kind of a typical commercial or hospital environment.

NOTE U_T is the AC mains voltage prior to application of the test level.

Guide & Manufacturer declaration - Electromagnetic immunity			
The device UDS-K\UDS-K LED is intended for use in the electromagnetic environment specified below. The customer or the user of the device UDS-K\UDS-K LED should make sure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guide
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3V 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the device UDS-K\UDS-K LED, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Minimum recommended separation distance $= [3,5/V_1] \times P^{1/2}$ $d = 1,2 \times P^{1/2}$ 80 MHz a 800 MHz $d = 2,3 \times P^{1/2}$ 800 MHz a 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol:
NOTE 1 From 80 MHz to 800 MHz the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device UDS-K\UDS-K LED is used exceeds the applicable RF compliance level above, the device UDS-K\UDS-K LED should be monitored to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device UDS-K\UDS-K LED. ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.			

**Recommended separation distances between
portable and mobile RF communications equipment and the device UDS-K\UDS-K LED**

The device UDS-K\UDS-K LED is intended for use in electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device UDS-K\UDS-K LED can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device UDS-K\UDS-K LED as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80MHz $d=1,2 \times P^{1/2}$	80MHz to 800MHz $d=1,2 \times P^{1/2}$	800MHz to 2,5GHz $d=2,3 \times P^{1/2}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) accordable to the transmitter manufacturer.

NOTE 1 From 80 MHz to 800 MHz the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

The product conforms to the following standards:

EN 60601-1:2006	EN 60601-1-2:2007
EN 60601-1-6:2010	EN 62366:2008
EN 61205:1994	EN ISO 22374:2005
EN 62304:2006	EN 980:2008
EN ISO 9687:1995	EN 1041:2008
ISO 15223-1:2012	EN ISO 14971:2012
EN ISO 7405:2008 +A1:2003	EN ISO 17664:2004
EN ISO 17665-1:2006	EN ISO 10993-1:2009
EN ISO 10993-5:2009	EN ISO 10993-10:2010

The device has been tested and homologated in accordance with EN 60601-1-2 for EMC. This does not guarantee in any way that this device will not be affected by electromagnetic interference. Avoid using the device in high electromagnetic environment.

11 – STATEMENT

All product modification rights are reserved to the manufacturer without prior notice. The images shown in this manual are for reference only.

The final interpretation rights belong to GUILIN WOODPECKER MEDICAL INSTRUMENT CO., LTD.

Design, interior parts, accessories etc. are patented and registered by WOODPECKER.

The production and sale of copies or fake products will be persecuted.

**TABLE OF OPERATING POWER OF THE TIPS**

	UDS-E (LED) UDS-L (LED) UDS-A (LED) UDS-P (LED)	UDS-J2 (LED)	UDS-J UDS-B UDS-K (LED)	UDS-N3 (LED)	UDS-N1 UDS-N2 (LED) UDS-N4	Compatible scaler brand
Scaling						
G1	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	Compatible with Woodpecker and EMS Ultrasonic Scalers
G2	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	
G3	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	
G4	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	
G5	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	
G6	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	
G7	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	
G8	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	
G9	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	
G10	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	
G11	1-10 (G)	1-10	1-9	LOW-HIGH	LOW-HIGH	
Periodontics						
P1	1-10 (P)	1-7	1-6	LOW-MID	LOW-MID	Compatible with Woodpecker and EMS Ultrasonic Scalers
P2L	1-3(P)	1-2	1-2	LOW	LOW	
P2LD	1-2(P)	1	1	LOW	LOW	
P2R	1-3(P)	1-2	1-2	LOW	LOW	
P2RD	1-2(P)	1	1	LOW	LOW	
P3	1-6(P)	1-4	1-3	LOW-MID	LOW-MID	
P3D	1-6(P)	1-4	1-3	LOW-MID	LOW-MID	
P4	1-6(P)	1-4	1-3	LOW-MID	LOW-MID	
Endodontics						
E1	1-3(E)	-	-	LOW	-	Compatible with Woodpecker and EMS Ultrasonic Scalers
E2	1-3(E)	-	-	LOW	-	
E3	1-6(E)	-	-	LOW	-	
E3D	1-3(E)	-	-	LOW	-	
E4	1-6(E)	-	-	LOW	-	
E4D	1-3(E)	-	-	LOW	-	
E5	1-6(E)	-	-	LOW	-	
E5D	1-3(E)	-	-	LOW	-	
P4D	1-6(E)	-	-	LOW	-	
E8	1-10(E)	-	-	LOW	-	
E9	1-10(E)	-	-	LOW	-	
E10	1-6(E)	-	-	LOW	-	
E10D	1-6(E)	-	-	LOW	-	
E11	1-6(E)	-	-	LOW	-	
E11D	1-6(E)	-	-	LOW	-	
E14	1-3(E)	-	-	LOW	-	
E15	1-3(E)	-	-	LOW	-	
Cavity Preparation						
SB1	1-10 (P)	1-7	1-6	LOW-MID	LOW-MID	Compatible with Woodpecker and EMS Ultrasonic Scalers
SB2	1-10 (P)	1-7	1-6	LOW-MID	LOW-MID	
SB3	1-10 (P)	1-7	1-6	LOW-MID	LOW-MID	
SBL	1-10 (P)	1-7	1-6	LOW-MID	LOW-MID	
SBR	1-10 (P)	1-7	1-6	LOW-MID	LOW-MID	

Note: "G" stands for Scaling Mode; "P" stands for Periodontics Mode; "E" stands for Endodontics Mode; "-" means not suitable for this model.

carlo de giorgi s.r.l.

Via Tonale, 1 – 20021 Baranzate (MI) – Italy
Tel. +39.02.356.15.43 r.a. – Fax +39.02.356.18.08
<http://www.degiorgi.it> E-Mail: info@degiorgi.it